

Shape Your Waist


# Shaken Massager

PEMF Smart Fat-Burning Device




A close-up photograph of a person's midsection. The person is wearing a black t-shirt and blue jeans. Their hands are pinching a portion of their belly fat. The background is plain white.

Is it challenging  
to lose fat?



It's **easy** when you  
**find the right**  
**method.**

The image features a blue wireframe cube on the left side, which is semi-transparent and appears to be floating. On the right side, there is a large, circular teal-colored graphic that resembles a complex data visualization or a stylized tunnel, composed of many concentric, slightly irregular lines and points, creating a sense of depth and movement. The background is a gradient of light teal to white.

To address  
the uncertain  
distribution of fat  
**through**  
**established physical laws**

# Smart Technology Leaves Fat Nowhere to Hide

## Shaken Massager PEMF Smart Fat-Burning Device



1. Ultrasound

Fat Reduction

2. Vibration

Body Shaping

3. Heat Therapy

4. Massage


5. RF  
(Radio Frequency)

Waist and Abdomen Toning

6. Red Light

7. PEMF

## Fat Reduction: Intelligent Sonic Technology to Target Stubborn Fat



### Ultrasonic

Utilizing CellBreak™ Cell Disruption Technology, it precisely targets the subcutaneous lipid layer (10-20mm deep), employs ultra-high frequency resonance ( $\geq 40,000\text{Hz}$ ) to precisely induce resonant dissociation of adipocyte membranes, and catalyzes the breakdown of



### Vibration



High-frequency (TBD) vibrations stimulate the contraction of waist and abdominal muscles, compressing fat cells.

### Heat Therapy



Thermal penetration awakens abdominal tissues and cells, promotes blood circulation, and accelerates fat metabolism.


## Body Shaping: Build Waist and Abdominal Muscles, Say Goodbye to Sagging Fat

### Massage




Multiple airbags apply and release pressure, simulating kneading and massage motions by human hands, making waist and abdominal muscles firmer and contours smoother.

### RF (Radio Frequency)



High-frequency electromagnetic waves penetrate 3–5mm beneath the skin, generating heat to stimulate collagen regeneration and enhance muscle elasticity in the waist and abdomen.

### Red Light




Red light penetrates the epidermis, accelerates waste metabolism, reduces pigmentation and wrinkles, and helps improve stretch marks and surgical scars.

## Strengthen Your Core: Technology That Goes Beyond Fat Burning to Provide a “SPA” for Your Organs



**7.8Hz ultra-low-frequency PEMF cross-medium penetration technology delivers energy intervention up to 20cm deep into biological tissue. By regulating cell membrane potential and mitochondrial metabolic pathways, it significantly downregulates the expression of pro-inflammatory factors. It targets and alleviates various types of inflammation, such as enteritis, nephritis, prostatitis, gynecological inflammation, lumbar muscle strain, and lumbar disc herniation.**



**Deep Activation of Microcirculation and Cellular Vitality: Efficiently clears metabolic waste and simultaneously alleviates inflammatory responses.**

**Intelligent Enhancement of Intestinal Rhythmic Movement: Optimizes nutrient conversion efficiency that helps restore a natural and comfortable state, and improves various intestinal issues.**

## One-Button Switch for Three Energy Fields: Unlock the Ultimate Solution for Core Health



### Body Shaping Mode

#### 7-in-1 Targeted Fat Burning

(PEMF, Heat Therapy, Ultrasound, Massage, Vibration, RF Radio Frequency, Red Light)

**Function: Fat burning and body shaping, fitness support**



### Massage Mode

#### 5D Deep Repair

(PEMF, Massage, Vibration, Heat Therapy, RF Radio Frequency)

**Function: Strengthens the core and provides organ care**



### Heat Therapy Mode


#### 2-Effect Gentle Care

(PEMF, Heat Therapy)

**Function: Keeps the waist warm while supporting menstrual care.**

## Wrist-Controlled Precision: See the Results of Every Fat-Burning Session

Shaken Massager features intuitive buttons and a minimalist design. You can switch modes and adjust intensity easily, whether using the main unit's control panel or your smartwatch.



Easily access all functions with just two intuitive buttons.



### Smartwatch Control Interface


**Built-in Bluetooth:** Connects quickly via Bluetooth, supporting real-time operation in dynamic scenarios like running or yoga.

**Comprehensive Health Tracking:** Built-in sensors monitor key health metrics, including heart rate, blood oxygen and sleep quality.

## Wear Freely, Burn Fat with Elegance


### Dual-Screen Control at Your Fingertips

Real-time display of operational status: the left screen shows gear and intensity, while the right screen displays mode and countdown timer.




### Seamless Fit to Your Waistline

Featuring a joint with 70°-120° degrees of freedom for adjustment, it accommodates waistlines from 60 to 140 cm and is suitable for individuals weighing from 40 to 130 kg, ensuring it does not slip when bending or turning.



### Technology That Balances Strength and Comfort

Constructed from soft silicone material and precision cast 304 stainless steel powder, it exhibits stable properties, good thermal conductivity, wear and corrosion resistance, is hypoallergenic, and can be easily cleaned with an alcohol wipe.



# Suitable for



Sedentary individuals (to relieve lower back muscle strain)



Fitness and body-shaping enthusiasts



Individuals with abdominal obesity



Postpartum women (to improve the appearance of stretch marks)



Individuals with intestinal disorders



Those experiencing dysmenorrhea, irregular menstruation, or gynecological inflammation



Individuals with prostatitis, urinary urgency, or incomplete urination



Those with nephritis, cystitis, or similar conditions

# Measure Your Waist



- \*Waist Measurement Method:**
1. Use a soft tape measure with 0.1cm precision.
  2. Wear lightweight clothing, stand upright, look forward, and relax your abdomen.
  3. Wrap the tape measure around the thinnest part of your waist (1-2cm above the navel) or directly around the navel.
  4. Ensure the tape measure is snug against the skin, with enough room to fit two fingers comfortably.
  5. Record the measurement to the nearest 0.1cm.

World Health Organization (WHO) Standards for Abdominal Obesity:

**International Standard (applicable to Europe, Africa, and other regions):**

Category	Male	Female
Normal	<94cm	<80cm
<b>Abdominal Obesity</b>	<b>≥94cm</b>	<b>≥80cm</b>
High Risk	≥102cm	≥88cm

**Asian Population Standard (applicable to East Asia, Southeast Asia)**

分类	男性	女性
Normal	<90cm	<80cm
<b>Abdominal Obesity</b>	<b>≥90cm</b>	<b>≥80cm</b>
High Risk	≥95cm	≥90cm

## P90 Empowered, More Efficient Fat-Burning Through Synergy



### **P90:**

Static aerobic exercise  
Activates cellular energy stations  
Burns full-body fat




### **Shaken Massager:**

Targets and breaks down abdominal fat  
Fat-burning and body shaping  
Improves abdominal and waist health

Twice a day, with at least 4 hours between sessions. Do not use simultaneously.  
It is recommended to use Shaken Massager first.

## Tech-Powered Fat-Burning, Easy and Comfortable



For professional use



For recreational use



For athletic activities



For business trips



# Shaken Massager

PEMF Smart Fat-Burning Device

Shake off fat and shine with confidence every day!